



2525 Davis Blvd, Unit B • Naples, FL 34104 • (239) 919-8719 • rdcmillling@gmail.com

593 SW Biltmore St, Port Saint Lucie, FL 34983 • (772) 785-9987 • rdcmillling@gmail.com

Sending ITERO Digital Impressions

We support our current and potential customers with the ability to accept digital impressions from all digital impression systems manufacturers. Digital impressions have several advantages over conventional impressions. They provide you with the ability to view your preparations, scan full arch cases and to communicate more effectively with your dental laboratory.

Instructions:

ITERO

1. Create New Case
2. Fill out treatment information (patient info, scanning order, restoration type, and shade). Choose RDC labs as lab of choice. To set up RDC Labs in your system, call ITERO customer support. Our account number is 63403.
3. Scan Patient.
4. Review the scan quality, interocclusal clearance, and margins.
5. Click the “SEND” icon in the main toolbar to send the case to the Cadent Center.
6. To confirm that that the file has been sent (or is queued for sending), open the Case Manager.

If any issues arise, please contact ITERO Customer Service at (800) 577-8767. RDC Labs will be notified to approve your connection. Once this is completed, you will be able to send your cases to RDC Labs.

If you'd like us to walk you through the process, or if you have any questions, please don't hesitate to contact us at (239) 919-8791. Our CAD/CAM team will be happy to assist you!